



The Hall A Coordinate Detector at Jefferson Lab: Adventures in Construction, Assembly and Testing

Parker Reid

Saint Mary's University
Summer 2016




Jefferson Lab: A Brief Overview



- Located in Newport News, Virginia
- Particle Physics lab
- In the process of a 12GeV upgrade for its synchrotron accelerator
- Working specifically for Hall A (1 of 4 Halls)



Hall A and the Super Bigbite Spectrometer




Source: <https://physics.fiu.edu/research-groups/nuclear-physics-experiment/>

- Takes high energy electrons from the Continuous Electron Beam Accelerator Facility (CEBAF)
-  Bigbite is currently under construction
- Designed to function with the CEBAF upgrade to 12 GeV
-  Spectrometer to have a very high acceptance rate.


Electron Scattering in an Electron-Proton collision

- High Energy Electron collides with a Proton (Target) at some position
- Electron is emitted with a different energy and a scattering angle
- The goal is to determine the position of the collision and the scattering angle



The Coordinate Detector (CDET)

- Consists of 6 Modules
- 3 Modules on each side of mount (See figure)
- Each Module contains 28 scintillation bars
- All Scintillation bars angled towards expected collision area
- Expected vertical resolution of $\sim < 2\text{mm}$




CDET Continued

- Scintillation **bar** contains 14 thin scintillation paddles
- Fiber optic cable threaded through each paddle



Several Bundles of Fiber optic cables


- Every **bar** is connected to a Photomultiplier tube (**PMT**)
- The PMT outputs the signal that data acquisition requires
- ★ Problem! These PMTs are being **Re-used**. Some may not be fit for purpose



Thin Scintillation paddle

PhotoMultiplier Tubes

- Take advantage of high voltage and the photoelectric effect
- Single photons interact with the 1st layer and release electrons
- These electrons travel to the next layer and release more, etc. Until a usable signal is generated



Source: www.franksweb.space.org.uk/

What we noticed:


1. A single PMT only has 1 High Voltage connection
2. Each scintillation bar contains 14 paddles, while each PMT contains 16 usable pixels


This still gives an opportunity to optimize!

Front and back of PMT used in CDET




PMT Data: Summer 2015


- Using PMT data from Nathan Murtha during summer 2015
- For each pixel:  Measured output current given a constant voltage and photon source
- Mean currents compiled into a text document for analysis




Example of data for a single pixel



Optimizing the Photomultiplier Tube Output Signal

- 2 Pixels with the largest deviation from the mean were removed 
- The ratio of the largest and smallest current of the remaining pixels is calculated
- Define as Maximum Scaling Factor
- Use as a basis for optimal PMT selection
- Module maps followed

Maximum Scaling Factor of each PMT (sorted)



		TOP LEFT		Notes	Notes	TOP RIGHT		Notes
Pixels NOT Used	Bar Label	PMT SERIAL #				PMT SERIAL #	Bar Label	Pixels NOT Used
1, 15	M3-8DL	154				52	M3-8DR	4, 14
1, 4	M3-9DL	102				108	M3-9DR	2, 16
2, 6	M3-10DL	33				99	M3-10DR	1, 4
2, 6	M3-11DL	31				107	M3-11DR	13, 16
14, 16	M3-12DL	150				181	M3-12DR	1, 4
1, 13	M3-13DL	46		MODULE 3		73	M3-13DR	13, 16
1, 4	M3-14DL	35				38	M3-14DR	3, 16
1, 13	M3-15DL	65				159	M3-15DR	4, 14
4, 12	M3-16DL	152				1	M3-16DR	1, 14
6, 13	M3-17DL	51				12	M3-17DR	5, 14
11, 14	M3-18DL	53				39	M3-18DR	2, 3
2, 13	M3-19DL	131				42	M3-19DR	4, 15
7, 11	M3-20DL	72				129	M3-20DR	4, 16
2, 4	M3-21DL	83				165	M3-21DR	1, 10
		BOTTOM LEFT				BOTTOM RIGHT		
		PMT #				PMT #		




High Voltage System: CDET HV
CAEN SY527N




High Voltage Control

- Previously no dedicated HV for CDET
- Designed HV test cables
- New HV switchboard
- Designed a map from HV remote control GUI to CDET
- Control on Linux based computer, operable from anywhere




Construction of the Coordinate Detector




Scintillator installation

- ❖ Tested with cosmic rays prior to installation (top left)
- ❖ Rejected pixels selected from prior analysis (top right)
- ❖ Fiber optic cables threaded through each bar
- ❖ Bars installed on bottom panel (bottom)




Installing Plates

- ❖ Aligned scintillation bars (top right)
- ❖ Installed top and side plates
- ❖ Taped overlapping areas with Electric tape (bottom right)



Finishing Touches

- ❖ Connect HV cables to PMT's
- ❖ Install discriminator cards
- ❖ Use HV to probe for light tightness
 $700\text{V} \rightarrow <30\text{nA}$
- ❖ Connect the remaining signal cables to data acquisition modules



Thank you

Questions?
